Rev. 6.0

Deutsch

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1. Product information:

The K3Pro® Implant System is a self-contained system of endosseous implants with corresponding surgical and prosthetic parts and instruments. K3Pro® Implants are made of titanium grade 4. The cover screws of the K3Pro® Implant System are made of titanium grade 5. This material has been repeatedly found to be biocompatible in studies. The surface configuration has been carefully designed and should not be altered. The implants are supplied with the OsteoActive® etched surface treatment.

Group	Name
ImplantImplantsate	Rapid/R-Line Sure/S-Line Compress/C-Line SHORT Rapid BS
Basic UDI ++EAGNIMPLANTSPK	

1.1 Intended use:

The K3Pro implant is intended for placement in the bone as a submucosal, transcutaneous or transcutaneous screw to which a permanent or temporary prosthesis can be attached. The K3Pro dental implant is designed for use in edentulous sites in the mandible or maxilla to support complete dentures, a terminal or intermediate abutment for fixed bridges or for partial dentures or as a single tooth replacement.

K3Pro dental implants can be used in both one- and two-stage surgery and their use is at the discretion of the clinician.

Immediate loading of K3Pro dental implants is at the discretion of the clinician and is dependent on bone quantity and quality, with non-occlusal loading and splinting being an important factor.

Immediate loading of K3Pro dental implants is at the discretion of the dentist and depends on bone quantity and quality, with non-occlusal loading and splinting being an important factor.

1.2 Intended users:

Surgeons, dentists and medical doctors familiar with and trained in augmentation.

1.3 Patient population:

Patients with completed alveolar growth.

Products for common use

The K3Pro implants are designed to be used together with the K3Pro abutments, the K3Pro transfer accessories & the K3Pro instruments.

For more detailed product descriptions including possible combinations and corresponding article numbers, please refer to the product catalogs on our website: www.argon-medical.com

2. Application:

The K3Pro® implant system can be placed in the mandible or maxilla. There are different sizes to choose from, which fit into healed or fresh extraction sites. The K3Pro® implants must be placed 1-2mm subcrestally, except for the K3Pro® SHORT implants (see 8.). For information on surgical procedures, please visit our website at www.argon-medical.com or take a course on K3Pro® Implant System surgical procedures.

3. Indications:

The K3Pro® Implant System is intended for use in edentulous sites in the mandible or maxilla to support a complete denture, bar, intermediate abutment for fixed bridges, partial dentures or as a single tooth replacement. The K3Pro® Implant System is intended for placement in bone as a submucosal screw to which a permanent or temporary prosthesis can be attached.

4. Contraindications:

The K3Pro® Implant System should not be used in patients with contraindicated conditions such as blood dyscrasias, uncontrolled diabetes, hyperthyroidism, bruxism, oral infections or malignancies. The K3Pro® Implant System should also not be used in patients with contraindicated conditions such as myocardial infarction within the last year or insufficient supporting bone to allow the use of an appropriately sized implant. Implants should not be placed if the width and height of the alveolar bone is insufficient to completely surround the implants.

a) K3Pro® Short Implants are not intended to be used as single-tooth replacements, neither splinted nor unsplinted, in the anterior tooth area, premaxillary or molar region. They should also not be used as denture fixation with less than four (4) implants. Their use is as a fixation for bridges with a maximum of one pontic between the abutments, in conjunction with multiple implants.

b) K3Pro® Mini Implants (ø 3.0 - 3.5 mm) are not intended for single tooth replacement, either splinted or splintless, in the premaxillary or molar region. They are also intended as permanent support for a prosthesis, but only when four (4) implants or more are used. They are not designed for lateral loading. Their only use is as an abutment for fixed bridges in conjunction with multiple implants.

5. Warnings:

a) Implant surgery is a highly specialized and complex procedure. The established techniques of oral implantology require special training. Implant courses and seminars are strongly recommended before surgery is performed. Incorrect technique can lead to implant failure and significant loss of surrounding bone. The K3Pro® Implant System should not be used in locations or situations other than those specifically indicated.

This may lead to failure of the implant with simultaneous destruction of the supporting bone. The K3Pro® Implant System may only be used with the bone drills specifically supplied for this purpose. Only dentists and oral and maxillofacial surgeons or oral surgeons are permitted to use the K3Pro® Implant System products.

b) Hygiene and disinfection: K3Pro® implants are supplied sterile. If the implant package is damaged or if the implants become contaminated in any way, they should be discarded. The K3Pro® Implant System is intended for single use. Implants are delivered sterile by irradiation and should not be

cleaned or sterilized.

c) MRI Safety Information:

A patient with a K3Pro Implant System implant can be scanned safely in an MRI system under the following conditions: Static magnetic field of 1.5 Tesla and max. 3 Tesla. Failure to comply with these conditions may result in injury. Based on currently available scientific literature, a temperature increase of <1°C can be expected for titanium dental implants, taking into account the magnetic field strength mentioned. Minor susceptibility artefacts may occur.

d.) Reporting of incidents:

All serious incidents that may occur in connection with the Argon Medical Dental Implant System in a country of the European Union must be reported to the manufacturer and the competent authority in your country. A serious incident means an incident that directly or indirectly caused, could have caused or might have caused any of the following:

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration in the state of health of a patient, user or other person,
- a serious risk to public health.

This requirement applies equally to users and patients. Please inform the patient of this requirement. The contact information for your country can be found on the website of the European Commission.

6. Precautions:

Panoramic radiographs or X-ray should be used to determine if sufficient bone is available at the intended implant site and to locate established critical anatomic features such as the mandibular canal, mental foramen, maxillary sinuses, and adjacent teeth. Adequate reporting and direct visual inspection of the prospective implant site are also required to determine the anatomy of the available bone. A thorough clinical evaluation is essential. Proper patient motivation is necessary if the procedure is to be successful.

7. Side effects:

The following complications associated with the procedure may occur: dehiscence, delayed healing, paresthesias, edema, hematoma, infection, inflammation, and general allergic reactions.

8. Notes on the K3Pro® Short implant system:

When using K3Pro® Short implants, the implant should be placed at least 3 mm subcrestally. Clinicians should carefully monitor patients for any of the following conditions: peri-implant bone loss, changes in implant response to vibration, or radiographic changes in bone at implant contact along the implant length. If the short implant has mobility or more than 50% bone loss, the implant should be evaluated for possible removal. If the clinician chooses a short implant, clinicians should consider a two-stage surgical approach in which a short implant serves as an additional implant and is placed with the largest possible implant diameter. Allow extended periods for osseointegration and avoid immediate loading.

9. Sterile handling:

The assistant removes the sterile Tyvek blister pack from the box, carefully opens the Tyvek blister pack and allows the inner pack to fall freely onto a sterile tray. Do not contaminate the inner packaging.

The assistant or dentist grasps the K3Pro® implant container and removes the cover cap, which also contains the K3Pro® cover screw. The implant is removed from the container in the following way:

a) Implant insertion:

Grasp the container with sterile gloved fingers and insert the insertion tool into the implant while squeezing the container. The insertion tool should engage the hexagonal connection of the implant, which is also indicated on the instrument mark. If inserted correctly, the mark on the insertion tool should meet the upper part of the implant body. The implant should now be slightly tapered and allow placement into the osteotomy without having to separate the implant and instrument. The implant is now placed into the prepared osteotomy according to the surgical instructions from the K3Pro® Implant System.

b) Instrument removal:

Remove the insertion instrument from the implant by gently wiggling the insertion instrument until the conical connection disengages from the implant and the instrument can be removed.

c) Placement of the cover screw of the K3Pro® implants:

Remove the cover screw from the container cover cap using a 1.2mm Hex screwdriver and turn counterclockwise. The screwdriver has a conical tip so that the connection with the cover screw is ensured for placement. Used correctly, there is no direct contact with the sterile component.

10. Surgical Instructions:

The implant bed is prepared according to the implant insertion length with a 2.0 mm twist drill using external irrigation. The implant bed is prepared with the step drills (3.0 - 6.0) up to the implant diameter at low speeds (40 - 50 rpm for titanium drills) without irrigation, (600-1000 rpm for steel). After completion, the instrument is center-punched to match the corresponding implant diameter.

For D1 - D2 type bone densities, the use of a tap (20 rpm) is recommended. It should be ensured that the implant length has been properly reamed and is not obstructed by bone particles. Safety can be ensured by using the bone level gauge. The K3Pro® implant is screwed into the prepared alveolus manually or by machine.

11. Recommended torques for implants:

Platform	troques (ratchet)	Max. insertion speed (hand-piece/-contra-angle)
2 mm	35 N/cm2	20 U/min
3 mm/ Short	45 N/cm2	20 U/min

The implant is intended to be used as a support for a temporary or permanent prosthesis, so these should be closed with a closure of the mucosal flap during osseointegration of the endosseous implants after implant placement.

12. The K3Pro® implant system warranty and limitations:

The success of all dental implants depends on the surgical procedure, which must be performed carefully. Careful patient selection and final restorative prostheses that meet individual patient needs and remaining anatomical conditions are essential. Selection of single or multiple implants of appropriate size and configuration for the existing anatomy is also critical to success. For the above reasons, Argon strongly recommends that all dentists attend one or more endosseous dental implant courses or symposia on the use of the K3Pro® Implant System before attempting to place implants in their dental practices. Because Argon cannot control factors within the scope of services performed by the dentist, including patient selection, surgical technique and restorative techniques, Argon assumes no responsibility beyond product replacement for failures or other adverse reactions or results from its purchasers arising from or through the use of a K3Pro® Implant. The K3Pro® Implant System is manufactured from Grade 4 surgical titanium and is supplied with the OsteoActive® surface treatment and, which has been found to be highly biocompatible by biomaterials researchers in numerous animal and human studies conducted over many years.

The K3Pro® Implant System and titanium abutments are exceptionally resistant to bending and fracture. However, there is no guarantee that fracture or ending of an implant or abutment cannot or will not occur due to traumatic injury or excessive occlusal stress. For the foregoing reasons, unless otherwise expressly stated, Argon makes no warranties beyond those contained in this warranty or beyond the description contained in any invoice. This warranty is

in lieu of all other warranties, expressed or implied, including merchantability and fitness for a particular purpose. Argon sells and warrants this to purchasers who are licensed dentists and who purchase implants and abutments for the purpose of their use and resale as part of the services they provide to patients. Argon does not grant any written warranty to any consumer, patient or end user and does not authorize any person to make such a written warranty on its behalf. The remedies of any purchaser from Argon (and its liabilities) shall be limited solely to the replacement or repair of a defective K3Pro® implant or abutment if such product was defective at the time of shipment and was returned to Argon for inspection, replacement or repair within six (6) months of its shipment. Argon shall not be liable to any person for incidental, consequential or other damages of any kind, whether a claim or potential claims are based on contract, negligence or tort theories.

In accordance with the Medical Device Regulation (EU) 2017/745 (MDR), Annex I, section 23.4 d), we ensure that the Summary of Safety and Performance evaluation (SSCP) of the device is accessible to the user and/or patient on our website.

13. Residual risks:

Despite all safety precautions, certain residual risks may exist when using this product:

- Risk of infection
- Rejection reactions
- Healing complications
- Implant displacement or failure
- Long-term effects

These risks can be minimized by correct use and regular medical examinations.

14. Safe disposal of the product:

Although not designed for disposal after implantation, it must be considered in accordance with the applicable legal requirements and the requirements of the notified body, as in the case of removal or other circumstances where the product is no longer required in the body. If the implant is removed or is no longer required for other reasons, please observe the following instructions for safe disposal:

Do not dispose of with household waste: This product must not be disposed of with normal household waste. It is necessary to dispose of the implant after removal via appropriate disposal channels specifically designed for medical waste.

Disposal via authorized disposal facilities: Contact an authorized medical waste disposal facility or your local waste disposal authority to dispose of the implant in a safe and environmentally sound manner.

Aster removal by a medical specialist: Should the implant be removed, disposal is carried out by the treating specialist or medical professional in compliance with the applicable regulations.

Please follow the applicable local regulations for the disposal of medical products to avoid risks to the environment and health.

CE 2797

Products of Argon Medical Productions & Vertriebs Gesellschaft mbH & Co. KG with CE approval meet the requirements of the Medical Devices Directive 93/42 EEC.

Ronly Caution: According to US Federal Law, this product may only be sold directly to or on behalf of a trained medical professional.

serial number

manufacturer

to be used until

warning

Single use only

do not use if package is damaged

sterilization by irradiation

product code

date of manufacturing

read instructions for use

Conditionally suitable for MRI

do not sterilize again

keep dry